

Automatic Reference Fuel Blending System for Octane & Cetane Number Determination



Shanghai Sinpar Scientific Instrument Co.,Ltd

Professional Manufacturer of Octane Engines

www.sh-sinpar.com

Automatic Reference Fuel Blending System

SINPAR FTC-AD1 Automatic Reference Fuel Blending Unit equipped with a modern computerized automatic control system, is used for the automated preparation of octane and cetane reference & standardization fuels blends for octane & cetane number testing according to ASTM D2699, ASTM D2700, and ASTM D613 standard test methods.

How to Accurately and Efficiently Prepare a Reference Fuel Blend?

In the traditional volumetric blending procedures, human error is unavoidable, and it is difficult to achieve both accuracy and speed.

Now, SINPAR FTC-AD1 Automatic Reference Fuel Blending System has solved this problem for you! With the automatic gravimetric measurement control system, it makes fuel blending *more accurate, faster* and *without human error*, and saves you a lot of time.

Why Choose FTC-AD1 Automatic Reference Fuel Blending Unit?

As is known to all, the accuracy of the reference fuel blend directly affects the result of the cetane/octane rating. A device that can prepare the fuel blend *accurately* and *efficiently* is very important to ensure the accuracy of octane/cetane number determination.

SINPAR is a professional manufacturer of fuel rating equipment. FTC-AD1 was developed based on SINPAR's many years of experience in fuel rating.

Over the years, it has been used not only in our own testing centers, but also in many laboratories around the world.

Its good and stable working performance has been well received by users. It is increasingly becoming a *must-have* in octane/cetane LABS.



Automatic Reference Fuel Blending System

Features & Benefits

■ Standards

ASTM D2699, ASTM D2700, ASTM D2885, and ASTM D613 standard methods.

■ High Precision

An electronic balance with the high precision of $0.01g$ is used for gravimetric measurement to ensure the accuracy of blending.

■ Simple Operation

A built-in touch-screen computer panel with *easy-to-use* installed software simplifies the blending operation.

■ Accuracy and Efficiency

The automatic calculation program controls the whole blending process to *eliminate human error* and improve the efficiency of octane/cetane testing.

■ Safety and Real-Time Monitoring

The unit is equipped with electrical protection & temperature monitoring system to meet labs safety requirements.

Fuel temperature and level in the tanks are monitored in real time, with *low fuel level & overfilling* warnings.


■ Reliability and Durability

The whole system is installed in a *movable and stable* industrial cabinet, which can connect to the laboratory's ventilation system.

The structure design ensures its reliability, safety, durability and usability.

Key Points

- ✓ high accuracy
- ✓ easy-to-use
- ✓ automatic control
- ✓ automated documentation
- ✓ save time
- ✓ economic
- ✓ safety
- ✓ durability



Automatic Reference Fuel Blending System

Automatic Blending By Weight

The blending unit prepares the reference fuel blends by weight (using a precision electronic balance) based on the ASTM test methods. With the automatic calculation and control program, it delivers the accurate and reliable fuel blend value with the precision of 0.01 ON/CN.




Professional Blending Software

The software was developed by SINPAR specifically for fuel blending. It is simple, clear, easy to use and requires no training.

From the main screen can access octane or cetane "BLENDING" operation, view blending "RECORD", perform pumps "CALIBRATION" and necessary "SETUP".


Status indicators of balance, vessel, tanks and valves are displayed in real time.



Octane/Cetane Blending Operation

The operator just needs to select the blending mode, enter the target blending parameters, and press "START" to begin the blending. The automatic computational control system controls the whole blending process. An accurate blend is delivered to the vessel automatically.

Free Blending mode is provided for special fuels blending.



Technical Specifications

Applicable Standards	ASTM D2699, ASTM D2700, ASTM D613
Application	For Octane and Cetane Number Determination
Measurement Mode	Gravimetric Measurement by Precision Balance
Fuel Blending Mode	Automatic Fuel Blending System with Software
Operating Mode	Built-in Touch Screen Computer
Tank Quantity	Optional Quantity(3/4/5/6) of Tanks with 10 L. Each
Blend Accuracy	±0.01 ON/CN
Blending Speed	1~2 min/500 ml
External Refilling System	Automatic Refilling Internal Tanks from External Tanks under Nitrogen Pressure (<i>option</i>)
Monitoring System	Real-Time Monitoring Fuel Temperature and Level; Low Fuel Level & Overfilling Warnings
Calibration Function	Precision Electronic Balance and Dosing Pumps
Blend Data	Auto-Saved Blending Certificate
Power Supply	100~240VAC 50Hz/60Hz with Single Phase
Safety System	Electrical Protection System & Temperature Monitoring
Ventilation System	Available to Connect with Aspirator or Ventilation Duct
Weight	200.00 kg (with six empty tanks)
Dimension	85.0x65.0x168.0 cm

**Due to continuing products development, the illustrations used may differ from actual products, and specifications are subject to change.*



Shanghai Sinpar Scientific Instrument Co., Ltd
Professional Manufacturer of Octane Engines

No.10, Yongchun East Road, Shanghai, China
T: +86-21-60538164 F: +86-21-60538100
sales@sh-sinpar.com www.sh-sinpar.com

©2022 SINPAR Company
All Rights Reserved